ACCESSORI VASCHE

TUB ACCESSORIES

GOC001

GOCCIOLATOIO SOLO SCHIENA

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SPLASH PROTECTION ONLY BACK

It allows to obtain maximum hygiene protecting wall from splashes and hairs during the washing, at the same time can be flanked on the sides by shelf for drying and mixer. Made of Aisi 304 stainless steel with S/B finish.

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The undersigned Lory Progetti Veterinarians srl, established in Via B. Buozzi 3/A, Cavriago, RE DECLARE That the equipment of which follows the characteristics: Denomination: GOC001

THIS ARTICLE DOES NOT REQUIRE INSTRUCTIONS BOOKLET

The appliance is built with materials corresponding to the European regulations for medical use in the veterinary field. The instrument does not contain any electrical parts and has a manual operation. Complies with the following directives: 2001/95/CE (general product safety) 2004/108/EC (Electromagnetic compatibility).

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TREATMENT INSTRUCTIONS FOR AISI 304 STAINLESS STEEL

FIRST CLEANING

Before using equipment made of aisi 304 stainless steel, it's recommended to clean the entire internal and external surface to remove any possible working residual; according to the type of equipment,

FOLLOW THIS INDICATIONS

- On a surface that will come into direct contact with water you must use a no-aggressive or anti-calcaree product. (Es like vialcal)
- On a surface that will not come everytime into direct contact with water you must use neutral

DETERGENT

Cleaning should be done by rubbing with a sponge and then rinsing thoroughly. For satin surfaces, clean following the satin finish. These instructions are necessary to avoid yellow spots that may surface due to a failure and correct cleaning before use.

PROTECTIVE PLASTIC

It is recommended, immediately after purchase, to remove the adhesive plastic from the surfaces that protects them, because, if leaved for a long time, it may release the residuals of glue which, not being visible immediately, may could turn yellow with the dirt deposit.

PERIODIC MAINTENANCE

For a good preservation of the stainless steel surfaces, especially when exposed to the atmosphere, it is sometimes advisable to resort to a periodical maintenance that consists of a washing with simple water and soap or neutral detergent taking care to provide to a good rinsing and to a subsequent careful drying. In any case it is always advisable to clean, rinse and dry the steel carefully with every use.

CONTACT WITH OTHER SURFACES

If other items are used above the steel surfaces, such as anti slip rubber mats, it is advisable to make sure that before placing them on the stainless steel surface, both surfaces are equally clean and dry with care. If this foresight is not implemented, in time the yellowing of stainless steel will be risked. At the end of the use, it is advisable to take care to remove the articles from the stainless steel surface when they are no longer used. The use of compounds that can release ionicloro, such as hydrochloric acid or its solutions, must be avoided, as stainless steel would be damaged.

THINGS TO AVOID:

- Use of wool or steel wool - vinegar vapours - hydrocyanic acid 100% - linoleic acid 100% up to 100 ° - Muriatic acid - sulphuric acid 5% boiling - sulphuric acid smoking oleum 50 ° - sulphurous acid 100% - Molten aluminium - sulphur dioxide 90% - sodium bisulfate 15% at 85 ° - saturated bleach - ferric chloride 5-50% - Ferrous chloride 10-20% - moist chlorine gas - calcium hypochlorite - sodium hypochlorite - Carbon tetrachloride - zinc chloride-ammonium phosphate - sodium phosphate - calcium hydroxide, Magnesium and potassium - borate sodium - sodium peroxide.