



P6

Waterproof Veterinary Ultrasound Scanner

Applications:

Ovine, Bovine, Equine, Swine, Alpaca, Feline, Canine, Rabbit, Fish, Snake and etc.

Features:

5.7" TFT-LCD, wide visual, high brightness, high contrast, image clear and exquisite.

Hand-held design, smart and light, easy to carry, convenient for detection.

Adopted AC and DC to supply power, lithium battery can be support machine to work more than 5 hours .

Specifications:


- ◆ Detect depth(mm): ≥ 140
- ◆ Blind zone(mm): ≤ 3
- ◆ Monitor size: 5.7 inch TFT-LCD
- ◆ Display mode: B, B+B
- ◆ Image gray scale: 256
- ◆ Image storage: available
- ◆ Cine loop: available
- ◆ Scan depth: available
- ◆ Image flip: up/down, left/right
- ◆ Focus position: adjustable
- ◆ Posture mark: available
- ◆ Image process: histogram, color encode, gama, image smooth
- ◆ Probe frequency: adjustable (3 point each probe)
- ◆ Frame correction: adjustable
- ◆ Measurement: distance, circumference, area, volume, GA and etc.
- ◆ Notation: date, time, name, sex, age, height, weight, farm.
- ◆ Output report: regular
- ◆ Charging method: adapter, battery
- ◆ Adapter rating: 100-240V, 50-60Hz, 70VA
- ◆ Battery continuous work: ≥ 5 hours
- ◆ Machine dimension: 203mm(L)x136mm(W)x45mm(H)
- ◆ Net weight: 950g

Standard Configuration:

- ◆ Scanner Body: 1pc
- ◆ Li-ion Battery: 1pc
- ◆ 3.5MHz Convex Probe: 1pc
- ◆ Charger with a Set of Wires: 1set
- ◆ Memory card: 8GB
- ◆ User Manual: 1pc

Optional:

- ◆ Micro-convex Probe
- ◆ Rectal Probe
- ◆ Convex Rectal Probe
- ◆ Back-fat Probe
- ◆ Li-ion Battery



Xuzhou Palmary Electronics Co., Ltd



ADD.: No.4 Building, Zhongkai City Light, West Huaihai Road, Xuzhou City, Jiangsu Province, China (221000)
 TEL.: 86-516- 83892858 / 83892898 FAX: 86-516-83892808
 BUSINESS: ultrasoundscanner@palmary.cn SERVICE: paierdianzi@palmary.cn
 www.palmaryultrasound.com