

# Surgery Covers

Surgical draping is the procedure of covering a patient and the surrounding areas with a sterile barrier to create and maintain a sterile field during a surgical procedure. The purpose is to eliminate passage of microorganisms between non-sterile and sterile areas.

## BUSTER Absorbent Surgery Cover

- Disposable 2-ply surgery cover
- Top layer of absorbent paper applied to polythene cover
- Absorbs fluid and blood, preventing fluids from leaking through
- Adhesive spray can be used for application



Cat. No	Description	Size	Package
141870	BUSTER Absorbent Surgery Cover	38 x 45 cm	25/pk
141871	BUSTER Absorbent Surgery Cover	45 x 75 cm	25/pk
141872	BUSTER Absorbent Surgery Cover	75 x 90 cm	25/pk

## BUSTER Surgery Cover

- Disposable surgery cover
- Flexible plastic
- 100% polyethylene
- Adhesive spray can be used for application



Cat. No	Description	Size	Package
141765	BUSTER Surgery Cover	30 x 45 cm	25/pk
141770	BUSTER Surgery Cover	60 x 90 cm	25/pk
141775	BUSTER Surgery Cover	45 x 60 cm	25/pk
141780	BUSTER Surgery Cover	90 x 120 cm	25/pk
141840	BUSTER Surgery Cover	120 x 120 cm	25/pk
141845	BUSTER Surgery Cover	120 x 250 cm	25/pk
141850	BUSTER Surgery Cover	120 x 180 cm	25/pk

## BUSTER Surgery Cover with adhesive edge around hole

- Absorbent surgery cover with adhesive edge around ellipse shaped hole
- Suitable for sterilisation and castration
- Non-woven top side and polyethylene film on the back side for fast liquid distribution and absorption
- The polyethylene film ensures that the drapes are impermeable to liquid and bacteria

### Principles of draping

**Isolation:** Isolates dirty, contaminated areas from clean areas

**Barrier:** Prevents fluid penetration

**Sterile field:** Creates a sterile surgical field by aseptic application of sterile materials

**Sterile surface:** Creates a sterile surface on the skin which acts as a barrier to help prevent skin flora from migrating to the incision site

**Fluid control:** Channels and collects body and irrigation fluids

**Equipment draping:** Allows non-sterile equipment to be brought into the sterile field



Cat. No	Description	Size	Package
141790	BUSTER Surgery Cover with adhesive edge around hole	XS	25/pk
141791	BUSTER Surgery Cover with adhesive edge around hole	S	25/pk
141792	BUSTER Surgery Cover with adhesive edge around hole	M	25/pk
141793	BUSTER Surgery Cover with adhesive edge around hole	L	25/pk
141794	BUSTER Surgery Cover with adhesive edge around hole	XL	25/pk
141795	BUSTER Surgery Cover with adhesive edge around hole	XXL	25/pk

## KRUUSE Adhesive Spray

- Adhesive spray used for application of BUSTER sterile covers
- Non-irritant and environmentally friendly, no C.F.C.



Cat. No	Description	Size	Package
141865	KRUUSE Adhesive Spray for Sterile Cover	150 ml	1/pk