




## VET LIFT

### Veterinary Surgery and Treatment Table

- High-end scissor type surgery and treatment table with versatile use in veterinary medicine
- Frame of maximum stability; powder coated according to RAL 7035
- Electrical, battery-powered or hydraulic height adjustment by foot switch from 31-101 cm
- Plate (stainless steel or synthetic top) with edge and drain (130 x 60 cm) as well as 14 tie-down hooks
- Optionally fixed with 2 castors installed on one side (allow smooth movement of the table by lifting on the opposite side) or moveable by 4 double castors with blocking mechanism
- Top optionally fixed or tiltable (continuously adjustable tilt angle of the top (up to 11°) with comfortable one-hand operation supported by gas spring)



#### Special version with integrated scale

- 0-100 kg in 50 g steps
- Animal weighing function (averaging method)
- Battery-operated (220 hours working hours)
- Automatic shutdown after 3 minutes
- Delivery without batteries (6 pieces AA)



Special version with integrated scale

## VET LIFT

fixed with small castors on one side



### Hydraulic height adjustable

Top	Item Code
fixed stainless steel top	769.605
tilting stainless steel top	769.606
fixed synthetic top	769.615
tilting synthetic top	769.616

### Electrical height adjustable (power cord)

Top	Item Code
fixed stainless steel top	759.605
tilting stainless steel top	759.606
fixed synthetic top	759.615
tilting synthetic top	759.616

### Battery-powered height adjustable

Top	Item Code
fixed stainless steel top	779.605
tilting stainless steel top	779.606
fixed synthetic top	779.615
tilting synthetic top	779.616

## VET LIFT

movable by 4 double castors (with blocking mechanism)



### Hydraulic height adjustable

Top	Item Code
fixed stainless steel top	769.601
tilting stainless steel top	769.602
fixed synthetic top	769.611
tilting synthetic top	769.612

### Electrical height adjustable (power cord)

Top	Item Code
fixed stainless steel top	759.601
tilting stainless steel top	759.602
fixed synthetic top	759.611
tilting synthetic top	759.612

### Battery-powered height adjustable

Top	Item Code
fixed stainless steel top	779.601
tilting stainless steel top	779.602
fixed synthetic top	779.611
tilting synthetic top	779.612

## VET LIFT with scale

fixed with small castors on one side



### Electrical height adjustable (power cord)

Top	Item Code
fixed stainless steel top	759.607
fixed synthetic top	759.617

### Battery-powered height adjustable

Top	Item Code
fixed stainless steel top	779.607
fixed synthetic top	779.617



Power supply: 230V 50/60Hz 1,5A

Lifting capacity: 120 kg