




BASE LIFT

Veterinary Surgery and Treatment Table

- Solid surgery and treatment table with versatile use in veterinary medicine made of highest-quality materials and components
- High-end stainless steel or synthetic top with edge and drain (130 x 60 cm)
- 14 tie-down hooks included
- Electrical (with power plug) or battery-powered height adjustment from 54 - 100 cm
- Top optionally stationary or tilting (continuously adjustable tilt angle up to approx. 30°) with comfortable one-hand operation supported by gas pressure spring
- Foot switch for height adjustment
- Optionally fixed with non-slip rubber feet or moveable by individually lockable double castors for easy maneuverability
- Comprehensive accessories available upon request
- Seamlessly pressed plate; completely lined underneath



BASE LIFT

fixed with non-slip rubber feet



Electrical height adjustable (power cord)

Top	Item code
tilting stainless steel top	638.600
tilting synthetic top	638.610
stationary stainless steel top	638.605
stationary synthetic top	638.615

Battery-powered height adjustable

Top	Item code
tilting stainless steel top	648.600
tilting synthetic top	648.610
stationary stainless steel top	648.605
stationary synthetic top	648.615

BASE LIFT

moveable by double castors (with blocking)



Electrical height adjustable (power cord)

Top	Item code
tilting stainless steel top	638.601
tilting synthetic top	638.611
stationary stainless steel top	638.606
stationary synthetic top	638.616

Battery-powered height adjustable

Top	Item code
tilting stainless steel top	648.601
tilting synthetic top	648.611
stationary stainless steel top	648.606
stationary synthetic top	648.616

Accessories for BASE LIFT

Article	Item code
OP rails (various)	602.596-602.600
Quick-release tie-down cleat*	602.601
Infusion pole*	602.602
Adjustable side supports (pair)*	602.604
Narcotic bow*	602.605
Leg support with tie-down cleat*	602.606
Tilting stainless steel dish with grid*	602.607
Rail extension (200 mm/pair)*	602.610
Instrument tray (300 x 200 mm)*	602.611
Arm rests for ophthalmology (pair)*	602.614
Top extension*	602.615, 602.616
Lateral positioners (pair)*	602.641, 602.646
X-ray translucent extension (90° tilting)*	612.232
Ultrasound table	602.640.1, 602.640.0-PB
Overlay mat black seamless (anti-bacterial)	70.123.053.02-P

*suitable for OP rails



Power supply: 230V 50/60Hz 3A

Lifting capacity: 120kg